ADEPT® Hip Resurfacing System

Patient Implant Information Leaflet

Read this leaflet carefully because it contains important information for you.

CONTENTS

ENGLISH (en)	3
FRANCAIS (fr)	14
DEUTSCH (de)	26
ESPAÑOL (es)	38
ITALIANO (it)	50
PORTUGUÊS (pt)	62
NEDERLANDS (nl)	74
ČESKÉ (cs)	86
MAGYAR (hu)	98
POLSKI (pl)	110
SLOVENŠČINA (sl)	122

ADEPT® Hip Resurfacing System

This leaflet provides information about:

- The type of medical device being considered;
- The type of medical condition the device is used for and how it works;
- The expected lifetime of the device and any necessary follow-up with the healthcare professional;
- Residual risks, precautions and measures to be taken if symptoms arise.

Device name

ADEPT® Hip Resurfacing System

Device model/type

ADEPT® / Hip implant

(Please refer to your implant card for the relevant product codes)

Intended purpose (use) of the device, indications for use and population in which the device is used

ADEPT® Hip Resurfacing System is intended to treat severely worn and diseased articulating surfaces of the hip joint caused by end-stage arthritis or trauma where non-surgical treatment routes have become less viable and the condition leads to pain, disability and compromised quality of life.

ADEPT® Hip Resurfacing System is intended for use in hip resurfacing arthroplasty for reduction or relief of pain and/or improved hip function due to osteoarthritis.

The ADEPT® Hip Resurfacing System is intended for skeletally mature patients.

Intended performance:

ADEPT® Hip Resurfacing System is intended to replace the worn or damaged cartilage within your hip joint. It is designed to work with the natural structures (ligaments and muscles) of your hip to provide you with a stable joint and a range of movement that allows you to perform normal everyday tasks with reduced or no pain. ADEPT® Hip Resurfacing System is expected to reduce pain and improve mobility compared to your pre-operative condition, as well as to have an acceptable risk of failure (revision) within your lifetime compared to other hip replacement devices.

Materials and substances included in the device

ADEPT® Hip Resurfacing System consists of two metal components (cobalt-chromium-molybdenum alloy). The femoral head is fixed to the thigh bone using a polymer cement (polymethylmethacrylate, PMMA) that cures when implanted. The acetabular cup (pelvis) component is press-fit into the bone which is shaped during surgery. A cobalt-chromium-molybdenum alloy beading and ceramic coating called hydroxyapatite on the acetabular cup promotes bone growth onto the implant.

There is overwhelming evidence that cobalt-containing implants do not represent an appreciable risk for systemic cobalt (or cobalt containing alloy) toxicity. The evidence, which is when exposed via clinically relevant routes, does not represent a hazard and supports that cobalt-containing alloys used in medical devices are unlikely to reach exposure levels that will elicit toxicity (polycythemia, thyroid changes, or myocardial damage), carcinogenicity, or reproductive toxicity. The favourable benefit-risk profile in relation to possible alternatives strongly support continued use of cobalt containing alloys in medical devices.

If you have any concerns about these statements, please speak to your doctor or healthcare professional.

Manufacturing residuals that could pose a risk to the patient

Manufacturing residuals are removed during validated cleaning processes before packaging and sterilisation of the product. No manufacturing residuals known to pose a risk to the patient remain.

Expected lifetime of the device

A well-performing device may last 10 years or more. However, as with any joint replacement device, individual or surgical factors can impact the device lifetime including, but not limited to:

- Female gender,
- Diagnoses other than osteoarthritis,
- Excessive activity,
- Failure to control body weight,
- Trauma such as a fall or impact to the joint,
- Abnormal loading or unreasonable stresses through your joint,
- Not following preoperative or postoperative rehabilitation advise,
- Significant postoperative leg length discrepancy (>1cm),
- Surgical technique including femoral neck notching, sub optimally sized implants, malalignment, inadequate fixation and/or cementing technique.

Special operating instructions for the use of the device, postoperative care (frequency of examination, monitoring, or maintenance of the device)

Operating instructions are provided to your healthcare provider and doctor for proper handling and implantation of the device. Talk to your doctor about self-care and physiotherapy in relation to your operation and the ADEPT® Hip Resurfacing System device.

There is no unique post operative care, monitoring or follow up required following an implantation of ADEPT® Hip Resurfacing System, however metal-on-metal bearings require blood tests for metal ion levels and/or additional imaging scans. You should be willing to attend regular follow-up clinics post-operatively and this may be required on an on-going basis.

Your doctor or healthcare provider will be able to advise you on postoperative care and physiotherapy in relation to hip resurfacing and this will be according to their normal practice of care.

Undesirable side effects and residual risks

The ADEPT® Hip Resurfacing System is not suitable for everyone and should not be used for the following reasons (contraindication to use):

- Active or suspected infection or sepsis in or around the hip joint or distant foci of infections (which may spread to the implant site).
- Osteoporosis, metabolic disorder, or a condition of marked bone loss which could impair bone formation.
- Cases where there is inadequate neuromuscular status, vascular insufficiency, muscular atrophy, incompetent or deficient soft tissue surrounding the joint, neuromuscular disease, poor skin coverage around hip joint or prior implantation that cannot provide adequate support or fixation for the prosthesis which would compromise implant stability or post-operative recovery, making the procedure unjustifiable.
- Skeletally immature patients.
- Severe dysplasia.
- Significant hip joint (acetabular socket or femoral bone) deformity that may prevent safe placement of components.
- Bone stock which in the surgeon's opinion is inadequate to support the device including bone absorption, severe osteopenia or a family history of osteoporosis, osteopenia or avascular necrosis. A Dual Energy X-Ray Absorptiometry (DEXA) scan may be necessary to assess inadequate bone stock.
- Number, position and size of cysts on the femoral bone head.
- Patients with known moderate to severe renal insufficiency.
- Patients who are immunosuppressed with diseases such as AIDS or persons receiving high doses of corticosteroids.
- Patients with known or suspected metal sensitivity (e.g. jewellery).
- Females of childbearing age or planning on becoming pregnant in the near term due to unknown effects on the foetus of metal ion release.

- Previous surgery or pre-existing metalwork that could compromise support, alignment or fixation and stability of the prosthesis which would compromise post-operative recovery, making the procedure unjustifiable.
- Patient has severe pre-existing conditions (comorbidities) or a limited life expectancy.

Residual risks / undesirable side effects include:

- Superficial or deep infection. Infections may occur months to years after surgery and these infections are difficult to treat and may require reoperation with revision surgery.
- Dislocation (separation of the articulating surfaces or partial separation known as subluxation) where muscle and fibrous tissue laxity (loose or weakened tissues around the joint) can contribute to these conditions.
- Accumulation of calcium around the joint or formation of bone (ossification) with loss of joint mobility, loosening or migration, wear or mechanical failure of the implant may result from trauma, loss of fixation, malalignment, bone resorption or deformation of the articular surface. Muscle and fibrous tissue laxity (loose or weakened tissues around the joint) or excessive loading can also contribute to these conditions.
- Bone remodelling (changes) are natural following joint replacement because the implant is a different material next to your own bone. Changes to the bone can lead to bone loss, implant loosening and failure.
- Loosening of the components can result in increased production of wear particles, as well as accelerate damage to the bone, making successful revision surgery more difficult.
- Poorly lubricated surfaces of the implant can lead to adverse wear and decreased function.
- Femoral neck fracture or bone perforation leading to bone fracture.
- Femoral head collapse and/or avascular necrosis (bone tissue dies due to lack of/interruption to blood supply).
- Fatigue fracture of an implant can occur as a result of excessive activity, malalignment, or trauma.
- (Heterotopic) Bone formation.
- Functionally reduced outcomes which may include increased pain and/or reduced hip function, temporary or permanent device related noise such as clicking or squeaking, leg length discrepancy, deformity.
- Delayed wound healing.
- Nerve pinching (impingement) or damage, vascular disorders (including thrombus) more frequently as a result of surgical trauma than the implant itself. Temporary or permanent nerve damage can result in functional and/or sensory deficits in the affected limb.
- Cardiovascular complications including venous thrombosis, pulmonary embolism or myocardial infarction and/or other systemic symptoms of cobalt toxicity.

- Gastrointestinal and/or urinary or genital complications.
- Sudden, pronounced, intraoperative drop in blood pressure due to the use of bone cement or other surgical factors.
- Hematoma (localised bleeding outside of blood vessel), or damage to blood vessels resulting in large blood loss intra-operatively.
- Adverse effects may necessitate reoperation, revision, arthrodesis of the involved joint, Girdlestone and/or amputation of the limb.

All of these adverse effects may require medical or surgical intervention or revision. Rarely these effects may lead to death.

Specific residual risks and undesirable side-effects related to metal-on-metal bearing are:

- Reactions due to metal particles may be described as (but not limited to) metallosis, adverse local tissue reaction (ALTR) and adverse reaction to metal debris (ARMD). Reactions to particles can cause bone resorption (osteolysis), other peri-prosthetic bone loss responses and loosening which may lead to the need to revise the components.
- Adverse reaction to metal debris/implant wear can result in histiocytosis granuloma formation, aseptic lymphocyte dominated vasculitis associated lesions (ALVAL), fluid collections, or soft tissue masses (pseudotumours), tissue damage and/or necrosis (tissue dies).
- Metal sensitivity or sensitisation including delayed Type IV sensitisation.
- Early and/or long-term increased serum, urine, and tissue levels of metal ions.
- Local reactions to metal debris can be asymptomatic. Regular follow-up is advised by the Health Authorities following any metal-on-metal procedure.
- Symptomatic patients may have pain, stiffness, joint instability, inflammation, and loosening around the hip joint. Local reactions should be identified early to avoid irreparable damage to the soft tissues.
- The risks of metal particles and ions, in particular cobalt, are not well quantified in the human body but studies have shown potential for cell damage, tissue necrosis, toxicity and inflammation. Local effects are more common but in rare cases systematic effects have been reported with a variety of symptoms often associated with hip pain. The presence of new or evolving hip symptoms in a patient warrants further investigation. Physicians and orthopaedic professionals should be aware of the wide range of signs and symptoms associated with systemic cobalt toxicity.
- Patients who have known metal sensitivity may be more prone to reaction to metal debris but currently there is no method to assess this before surgery.
- There is a risk that cobalt chromium alloys may be carcinogenic but there is no confirmatory scientific evidence for this.

Patients on medications (such as high dose of chronic amino glycoside treatment) or with pre-existing conditions (such as diabetes) that increase the risk of future, significant renal impairment should be aware of an increase in systemic metal ion concentration. Preoperative and postoperative monitoring of renal function may be necessary.

You should contact your healthcare professional if you believe that you are experiencing side-effects related to the device or its use, or if you are concerned about the risks. This document is not intended to replace a consultation with your doctor if needed.

While hip resurfacing devices are largely successful in attaining mobility and function, they do have a limited lifespan, and this is dependent on many factors that can be both surgical and patient related such as:

- Patients who are female.
- Component head size.
- Patients who are obese, failure to control body weight.
- Diagnosis other than osteoarthritis.
- Failure to follow post-operative care guidelines.
- Poor patients reported scores before surgery and/or at post-operative follow-up.
- Patients with significant leg length discrepancy (>1cm).
- Femoral neck (bone) notching during the surgical procedure.
- Loss of acetabular press-fit either during initial operation or post-operatively.
- Surgical factors such as suboptimal positioning, sizing or fixation of components and off-label use of components can lead to early failure by premature wear, abnormal loading or unusual stress, third body wear, or failure of the device. Your surgeon should be thoroughly familiar with the implants, instruments, and surgical procedure prior to performing surgery.

These risk factors are not exhaustive and current evidence should always be used by the surgeon when deciding if hip resurfacing is the right treatment choice for you.

The more risk factors a patient has, the greater the risk of the hip resurfacing device failing.

Levels of function and satisfaction will vary between patients and you should be aware the hip resurfacing may not restore function to the level expected with a normal healthy joint. Managing expectations before surgery can be helpful.

During the first post-operative 6 to 12 months the bone will be integrating with the device. You should discuss returning to leisure and sporting activities, which place higher demand on the joint replacement, with your surgeon who can offer the most up to date advice.

Patient guide for return to activity pathway after hip resurfacing include:

- Hip resurfacing patients should not return to activity faster than THR patients due to increased early postoperative risks, including femoral neck fracture, unexplained pain and early failure if activity is resumed too soon.
- Bone quality affects recovery – patients with lower bone density or related factors, may need extra precautions when returning to activity. Your surgeon will be happy to discuss your bone health and how this might affect your recovery and risks following surgery.
- All activities carry risks of accidental movements (stumbles, falls, excessive weight-bearing) that must be minimised in early healing.
- The bone healing process takes time and while resurfacing patients can achieve high activity levels long-term, patience in the early phase is crucial for success.
- 0–6 weeks (initial recovery & bone healing):
 - Weight-bearing: Allowed as tolerated but should be controlled. Use crutches or a walker initially.
 - Walking: Focus on smooth, balanced movement to avoid compensatory strain.
 - Restrictions: Avoid twisting, sudden weight shifts, or stumbles.
 - Daily Activities: Short, controlled walks and self-care tasks.
 - Risk Awareness: Unexplained pain is a sign to ease activity, not push through.
- 6 weeks – 3 months (gradual functional recovery):
 - Mobility & Strength: Transition to unaided walking but avoid excessive strain.
 - Physio Therapy: Focus: Strength training and balance exercises, with careful load progression.
 - Low-Impact Activities: Swimming, cycling, and elliptical training may begin.
 - Restrictions: Avoid twisting, sudden weight shifts, or stumbles.
 - Risk Awareness: New bone continues to form during this period, so increase activity cautiously. Unexplained pain is a signal to ease activity, not push through.
- 3–6 months (careful return to activity, no impact yet):
 - Exercise Expansion: Careful expansion of exercise such as controlled strength training and hiking on even terrain.
 - Restrictions: No running, jumping, or abrupt direction changes.
 - Risk Awareness: Continue to monitor how you feel and adjust activity if unexplained pain or discomfort arises. Bone is still forming during this period.
- 6–12 months (transition to full activity – with caution):
 - Higher-Demand Activities: e.g. Light jogging, skiing, or tennis may be considered with surgeon clearance.
 - Gradual Progression: Intensity and frequency should increase incrementally. Adjust activity if discomfort arises.
 - Mindful Movement: Avoid activities with a high risk of falls or sudden impact.
- Beyond 12 months (living a fully active life):
 - Most patients can resume a wide range of activities but should remain mindful of high-impact risks.

Warnings, precautions or measures with regard to exposure to reasonably foreseeable external environments including risks of device interaction with other equipment

The ADEPT® Hip Resurfacing has been tested in a laboratory for interaction with Magnetic Resonance Imaging (MRI) equipment. General recommendations for MR parameters which should not expose you to additional hazard during MRI are provided in the Instructions for Use which your healthcare provider can access here: '<https://info.matortho.com/reg/400-482.pdf>.'

You should make your doctor aware of any other implants or known health considerations that may require reduction of the MRI limits. There are several different manufacturers and generations of MRI systems available and MatOrtho® cannot make any claims regarding the safety of MatOrtho® implants and devices with any specific MRI system.

Symptoms that could indicate that the device is malfunctioning

Unexpected pain, stiffness, instability, or a general worsening of symptoms around your hip, especially when compared to your preoperative condition, may indicate that there is a problem with your joint replacement. An inflamed, hot, itchy, or red area around your hip or wound may indicate infection.

Precautions and measures that should be taken if the performance of the device changes or if any symptoms mentioned in this leaflet are experienced

A general worsening of the performance of the replaced hip joint as compared to the preoperative condition may indicate that there is a problem with the joint replacement. Serious incidents that occurs in relation to the device should be reported to the appropriate country authority:

- Australia: Therapeutic Goods Administration (TGA). Information for this can be found on the Therapeutic Goods Administration's website: <https://www.tga.gov.au/medical-device-incident-reporting-investigation-scheme-iris>.
- UK: <https://yellowcard.mhra.gov.uk/>
- EU (Authorised Representative):
 - > AR@matortho.com
 - > <https://www.hpra.ie/about-us/patient-forum/patient-forum-resources>
 - > <https://www.hpra.ie/safety-information/how-we-monitor-safety/medical-devices/reporting-safety-issues-for-devices/reporting-device-safety-issues-for-healthcare-professionals>

- Canada:
<https://www.canada.ca/en/health-canada/services/drugs-health-products/reports-publications/medeffect-canada.html>

The healthcare professional should inform the manufacturer of any serious adverse events by emailing 'per.admin@matortho.com.'

Precautions and other measures that should be taken at, or near, the end of the expected device lifetime

The precautions, risks, and warnings described above are appropriate throughout the lifetime of the device.

Other circumstances in which you should contact a health professional in relation to the operation of the device

No further specific circumstances beyond those discussed in this leaflet.

Manufacturer

For further information, please visit MatOrtho® patient information website:
www.matortho.com

or contact MatOrtho® at:

MatOrtho Limited, 19/20 Mole Business Park, Randalls Road, Leatherhead, Surrey, KT22 7BA, United Kingdom.

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