ReCerf™ Hip Resurfacing System

Patient Implant Information Leaflet

Read this leaflet carefully because it contains important information for you.

CONTENTS

ENGLISH (en)	3
FRANCAIS (fr)	10
DEUTSCH (de)	18
ESPAÑOL (es)	26
ITALIANO (it)	34
POLSKI (pl)	42

ReCerf® Hip Resurfacing System

This leaflet provides information about:

- The type of medical device being considered;
- The type of medical condition the device is used for and how it works;
- The expected lifetime of the device and any necessary follow-up with the healthcare professional;
- Residual risks, precautions and measures to be taken if symptoms arise.

Device name

ReCerf® Hip Resurfacing System

Device model/type

ReCerf® Hip Resurfacing/ Hip implant

(Please refer to your implant card for the relevant product codes)

Intended purpose (use) of the device, indications for use and population in which the device is used

ReCerf® Hip Resurfacing is intended to replace the articulating surfaces which make up the hip joint, for the treatment of pain and disability due to degeneration of cartilage at the bone interface and deterioration of the joint. Hip resurfacing is dictated by the disease progression and may be used when more conservative treatment routes have been exhausted.

ReCerf® Hip Resurfacing System is used to treat patients who have been diagnosed with arthritis of the hip. Treatment may also serve to correct a deformity or dysfunction of the hip resulting from the arthritis.

The intended purpose (or intended use) for ReCerf® Hip Resurfacing System is in hip resurfacing arthroplasty (replacement) for the reduction or relief of pain and/or improved hip function in skeletally mature adults having osteoarthritis.

The indication for use for ReCerf® Hip Resurfacing System is osteoarthritis of the hip.

Intended performance: ReCerf® Hip Resurfacing System is expected to reduce pain and improve mobility compared to the pre-operative condition, as well as to have an acceptable risk of failure (revision) within the patients' lifetime compared to other hip replacement devices.

Materials and substances included in the device

ReCerf® Hip Resurfacing System consists of two ceramic components made using BIOLOX® *delta* ceramic. One component replaces the surfaces of the femoral head and the other replaces the surfaces of the acetabular socket. Both these components articulate together to form the hip joint. The femoral head component is fixed using a polymer cement (polymethylmethacrylate, PMMA) and the acetabular cup component is press-fit (mechanical fixation) into the hip socket. The cup has a coating (titanium and hydroxyapatite) which helps to promote bone ingrowth directly onto the implant.

Manufacturing residuals that could pose a risk to the patient

Manufacturing residuals are removed during validated cleaning processes before packaging and sterilisation of the product. There are no known risks due to manufacturing residuals.

Expected lifetime of the device

ReCerf® is intended to have a lifespan of 10 years, which is in line with other hip resurfacing devices, however there is a lack of clinical evidence to support this. ReCerf® was first implanted clinically in September 2018 and the 2024 annual Australian Orthopaedic Association National Joint Registry (AOANJRR) published a revision rate of 0.6% (0.1, 2.3, 95% confidence intervals) at 3 years for ReCerf®. Survivorship beyond this timepoint cannot be predicted, however data and the length of time the device has been in clinical use is building and the most up to date information can be accessed on the AOANJRR website via <https://aoanjrr.sahmri.com/> and within the MatOrtho ReCerf® Clinical Brochure which can be found at www.matortho.com.

Risk factors discussed in this Patient Information Leaflet may limit the device lifetime. The more risk factors a patient has, the greater the risk of early failure of the hip resurfacing and subsequent revision of the prosthesis. These include pre-operative or progressive conditions that affect the quality of bone, excessive activity, unmanaged body weight, trauma such as a fall or impact to the joint, abnormal or uneven loading through the joint, failure to follow preoperative or postoperative rehabilitation advice, and surgical variables such as implant positioning, fixation, or anatomical changes like a significantly difference in leg length may also influence outcomes.

Special operating instructions for the use of the device, postoperative care (frequency of examination, monitoring, or maintenance of the device)

Your post-operative recovery is important and your doctor will be able to advise you on postoperative care and physiotherapy in relation to hip resurfacing and this will be according to their normal practice of care and there is no is no unique post operative care, monitoring or follow up required following a ReCerf® Hip Resurfacing device.

Operating instructions are provided to the healthcare professional for proper handling and implantation of the device.

Undesirable side effects and residual risks

While hip resurfacing devices are largely successful in attaining mobility and function, they do have a limited lifespan, and this is dependent on many factors that can be both surgical and patient related such as:

- Patients who are obese or excessive loading.
- Bone health.
- Intra-operative surgical factors such not following the recommended operative technique, femoral neck (bone) notching, incorrectly positioned or components not fully seated, loss of acetabular cup press-fit either during initial operation or post-operatively, poor cementing technique.

- Patients with significant leg length discrepancy (>1cm).
- A potentially higher risk of early failure for smaller component head sizes and females. This is a known risk factor for hip resurfacing devices made from metal materials.
- The more risk factors, the greater the risk of early failure of the hip resurfacing device and subsequent revision of the device.
- During the first post-operative 6 to 12 months the bone will be integrating with the device:
 - > Hip resurfacing patients should speak with their surgeon about their return to activity due to increased early postoperative risks, including femoral neck fracture, unexplained pain and early failure which can occur if activity is resumed too soon.
 - > Bone quality affects recovery and patients with lower bone density or related factors, may need extra precautions when returning to activity.
 - > All activities carry risks of accidental movements (stumbles, falls, excessive weight-bearing) that must be minimised in early healing.
 - > The bone healing process takes time and while resurfacing patients may achieve high activity levels long-term, patience in the early recovery phase is crucial for success.
- Infections may occur months to years after surgery and these infections are difficult to treat and may require reoperation with revision surgery.
- Dislocation (separation of the articulating surfaces or partial separation known as subluxation), loosening, migration, wear or mechanical failure of the implant, loss of fixation, malalignment, bone resorption, changes to the implant surfaces, accumulation of calcium or bone formation (ossification) or loss of joint mobility. rfaces, accumulation of calcium or bone formation (ossification) or loss of joint mobility.
- Trauma and fibrous tissue laxity (weak/loose tissue) or excessive loading can lead to poor outcome.
- Bone changes are natural following joint replacement because the implant is a different material next to the bone. Changes to the bone can lead to bone loss, implant loosening, bone fracture and failure.
- The coating on the acetabular shell component may separate from the ceramic implant component leading to loosening, pain, reduced functional outcome, revision.
- Loosening of components. Loosening can be caused by surgical factors such as component malalignment or suboptimal seating of components, poor cementation, increased wear due to abnormal function of the bearing surfaces and particles generated during wear may accelerated damage to the bone. Loosening can lead to pain, decreased function and early revision.
- Femoral neck fracture, notching of the bone, bone fracture, bone perforation leading to fracture, femoral head collapse and/or onset or progression of avascular necrosis (bone tissue dies due to lack of/interruption to blood supply). The pre-existing condition, reason for hip replacement, use of

medication/drugs or periods of limb disuse can affect bone quality.

- Fracture/breakage of an implant component may occur as a result of extreme or abnormal loading, malalignment of components, or trauma.
- (Heterotopic) Abnormal bone formation.
- Poor component positioning or suboptimal function can lead to increased wear, pinching (impingement), mechanical symptoms, unusual sounds, pain or reduced functional outcomes.
- Levels of function and satisfaction will vary between patients and the hip resurfacing device may not restore function to the level expected from a normal healthy joint. Managing expectations before surgery can be helpful. Functionally reduced outcomes may include increased pain and/or reduced hip function, temporary or permanent device related noise such as clicking or squeaking, leg length discrepancy or deformity.
- Delayed wound healing.
- Nerve pinching (impingement) or damage, vascular disorders (including thrombus) more frequently as a result of surgical trauma than the implant itself. Temporary or permanent nerve damage can result in functional and/or sensory deficits in the affected limb.
- General surgical risks which are not specific to the ReCerf® device include cardiovascular complications including venous thrombosis, pulmonary embolism or myocardial infarction; gastrointestinal and/or urinary or genital complications; sudden, pronounced, intraoperative drop in blood pressure due to the use of bone cement or other surgical factors; hematoma (localised bleeding outside of blood vessel), or damage to blood vessels resulting in large blood loss intra-operatively.
- Adverse effects may necessitate reoperation, revision, excision, arthrodesis of the involved joint, Girdlestone and/or amputation of the limb.

All of these adverse effects may require medical or surgical intervention or revision. Rarely these effects may lead to death.

Warnings, precautions or measures with regard to exposure to reasonably foreseeable external environments including risks of device interaction with other equipment

Magnetic Resonance Imaging (MRI) is a medical imaging tool that uses strong magnetic fields and radio waves to create detailed pictures inside the body.

The ReCerf® Hip Resurfacing device has not been tested for interaction with Magnetic Resonance Imaging equipment directly. However, a similar geometry metal-on-metal hip resurfacing device has been tested, and this represents a worse case than the ReCerf® device and the results inform the recommendations for Magnetic Resonance parameters. The patient is advised to discuss with the healthcare professional their ReCerf® implant and any other implants or known health considerations before undergoing Magnetic Resonance Imaging, because different factors may require reduction

of Magnetic Resonance Imaging limits. ReCerf® has a titanium metal coating and metal can cause artefacts during Magnetic Resonance Imaging. Metal artefacts can make it difficult to interpret the result of imaging. There are several different manufacturers and generations of Magnetic Resonance Imaging systems available and MatOrtho® cannot make any claims regarding the safety of MatOrtho® implants and devices with any specific Magnetic Resonance Imaging system.

Symptoms that could indicate that the device is malfunctioning

Unexpected pain, stiffness, instability, or a general worsening of symptoms around your hip, especially when compared to your preoperative condition, may indicate that there is a problem with your joint replacement. An inflamed, hot, itchy, or red area around your hip or wound may indicate infection.

Precautions and measures that should be taken if the performance of the device changes or if any symptoms mentioned in this leaflet are experienced

A general worsening of the performance of the replaced hip joint as compared to the preoperative condition may indicate that there is a problem with the joint replacement. Serious incidents that occurs in relation to the device should be reported to the appropriate country authority:

- Australia: Therapeutic Goods Administration (TGA). Information for this can be found on the Therapeutic Goods Administration's website: <https://www.tga.gov.au/medical-device-incident-reporting-investigation-scheme-iris>.
- UK: <https://yellowcard.mhra.gov.uk/>
- EU (Authorised Representative):
 - > AR@matortho.com
 - > <https://www.hpra.ie/about-us/patient-forum/patient-forum-resources>
 - > <https://www.hpra.ie/safety-information/how-we-monitor-safety/medical-devices/reporting-safety-issues-for-devices/reporting-device-safety-issues-for-healthcare-professionals>
- Canada:
<https://www.canada.ca/en/health-canada/services/drugs-health-products/reports-publications/medeffect-canada.html>

The healthcare professional should inform the manufacturer of any serious adverse events by emailing 'per.admin@matortho.com.'

Precautions and other measures that should be taken at, or near, the end of the expected device lifetime

The precautions, risks, and warnings described above are appropriate throughout the lifetime of the device.

Other circumstances in which you should contact a health professional in relation to the operation of the device

No further specific circumstances beyond those discussed in this leaflet.

Manufacturer

For further information, please visit MatOrtho® patient information website:
www.matortho.com

or contact MatOrtho® at:

MatOrtho Limited, 19/20 Mole Business Park, Randalls Road, Leatherhead, Surrey, KT22 7BA, United Kingdom.

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