# **SAIPH® Knee System**

## **Patient Implant Information Leaflet**

Read this leaflet carefully because it contains important information for you.

## CONTENTS

---

<b>ENGLISH</b> (en)	<b>3</b>
<b>FRANCAIS</b> (fr)	<b>11</b>
<b>DEUTSCH</b> (de)	<b>19</b>
<b>ESPAÑOL</b> (es)	<b>27</b>
<b>ITALIANO</b> (it)	<b>35</b>
<b>PORTUGUÊS</b> (pt)	<b>43</b>
<b>NEDERLANDS</b> (nl)	<b>51</b>
<b>ČESKÉ</b> (cs)	<b>59</b>
<b>MAGYAR</b> (hu)	<b>67</b>
<b>POLSKI</b> (pl)	<b>75</b>
<b>SLOVENŠČINA</b> (sl)	<b>83</b>

# SAIPH® Knee System

This leaflet provides information about:

- The type of medical device being considered;
- The type of medical condition the device is used for and how it works;
- The expected lifetime of the device and any necessary follow-up with the healthcare professional;
- Residual risks, precautions and measures to be taken if symptoms arise.

## Device name

SAIPH® Knee System

## Device model/type

SAIPH® Knee / Knee implant

*(Please refer to your implant card for the relevant product codes)*

## Intended purpose (use) of the device, indications for use and population in which the device is used

SAIPH® Knee is intended for use in total knee replacement for the replacement of the articulating surfaces of the knee joint (which includes the tibiofemoral and patellofemoral joint) for the treatment of pain and disability due to arthritis.

The indications (symptoms suggesting a medical treatment is required) for which SAIPH® Knee is used are osteoarthritis including post-traumatic osteoarthritis, rheumatoid arthritis and previously unsuccessful knee replacement, where one or more compartments of the knee are involved. The SAIPH® Knee System is intended for skeletally mature patients.

## Intended performance:

SAIPH® Knee System is intended to replace the worn or damaged cartilage within your knee joint. This includes the cartilage on the end of your thighbone (femur), the top of your shin bone (tibia) and, optionally, the cartilage on the underside of your kneecap (patella). SAIPH® Knee System is designed to work with the natural structures (ligaments and muscles) of your knee to provide you with a stable joint and a range of movement that allows you to perform normal everyday tasks with reduced or no pain.

## Materials and substances included in the device

The SAIPH® Knee System consists of two metal components (cobalt-chromium-molybdenum alloy) and polyethylene (ultra-high molecular weight polyethylene, UHMWPE) component(s). Both metal components are fixed to the bone (one fixed to the femoral bone and one fixed to the tibial bone of the knee joint) using a polymer cement (polymethylmethacrylate, PMMA). The UHMWPE component is fixed directly and securely to the tibial metal component. The UHMWPE component

then articulates with the metal femoral component during movement of the knee joint. The kneecap (patella) may be replaced with a UHMWPE component, at the discretion of your operating surgeon and when selected it is implanted also using PMMA cement. The SAIPH® Knee performs well with and without a replaced patella.

There is overwhelming evidence that cobalt-containing implants do not represent an appreciable risk for systemic cobalt (or cobalt containing alloy) toxicity. The evidence, which is when exposed via clinically relevant routes, does not represent a hazard and supports that cobalt-containing alloys used in medical devices are unlikely to reach exposure levels that will elicit toxicity (polycythemia, thyroid changes, or myocardial damage), carcinogenicity, or reproductive toxicity. The favourable benefit-risk profile in relation to possible alternatives strongly support continued use of cobalt containing alloys in medical devices.

If you have any concerns about these statements, please speak to your doctor or healthcare professional.

### **Manufacturing residuals that could pose a risk to the patient**

Manufacturing residuals are removed during validated cleaning processes before packaging and sterilisation of the product. No manufacturing residuals known to pose a risk to the patient remain.

### **Expected lifetime of the device**

A well-performing device may last 10 years or more. However, as with any joint replacement device, individual or surgical factors can impact the device lifetime including, but not limited to:

- Excessive activity,
- Failure to control body weight,
- Trauma such as a fall or impact to the joint,
- Abnormal loading or unreasonable stresses through your joint,
- Not following preoperative or postoperative rehabilitation advise,
- Poor surgical technique and/or deviation from the SAIPH® prescribed surgical technique including sub optimally sized implants, malalignment and/or soft tissue imbalance which can place unusual forces on the components, inadequate fixation and/or cementing technique, or inadequate wash-out prior to wound closure which may leave residual particles from cementing, bone fragment or other debris. Your surgeon is required to be thoroughly familiar with the implant, instruments, and surgical procedure prior to performing surgery.

### **Special operating instructions for the use of the device, postoperative care (frequency of examination, monitoring, or maintenance of the device)**

Operating instructions are provided to your healthcare provider and doctor for proper handling and implantation of the device.

There is no unique post operative care, monitoring or follow up required. However, your surgeon or healthcare provider will be able to advise you on postoperative care and physiotherapy in relation to knee replacements and this will be according to their normal practice of care.

## **Undesirable side effects and residual risks**

As with any knee replacement procedure residual risks and undesirable side effects can occur. This leaflet makes you aware of the surgical and device related risks and makes you aware of the possible adverse effects. These include:

- Prolonged stay in hospital.
- Superficial infection, cellulitis or deep infection; infection can lead to the failure of the joint replacement.
- Foreign materials and debris from the implant can lead to inflammation (histiocytosis granuloma formation) or adverse reaction such as sensitivity (which may lead to consequent osteolysis and loosening of the implant).
- Loosening or migration, due to trauma, loss of fixation, malalignment or bone resorption.
- Soft tissue imbalance and/or malalignment of the components which can cause excessive wear and/or failure of the implant.
- Dislocation (separation of the articulating surfaces or partial separation known as subluxation) due to improper positioning of the implant components which leads to postoperative instability. Muscle and fibrous tissue laxity (loose, weak tissue surrounding the joint) can also contribute to these conditions.
- Fatigue fracture of the implant due to excessive activity, malalignment, or trauma.
- Instability or implant subsidence (the implant may move within the bone around the implant).
- Wear or deformation of the articular surface as a result of excessive loading.
- Soft tissue damage or rupture(s).
- Stiffness, arthrofibrosis (abnormal growth of fibrous tissue around the joint) or suboptimal range of motion which may lead to further surgery such as manipulation under anaesthetic (MUA) or revision.
- Bone fracture during surgery or postoperatively due to trauma, for example fall, or excessive loading, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, or bone resorption.
- Bone notching, bone perforation or bone fracture.
- Patella (kneecap) impingement, additional operations on the patella including the addition of a patella if you did not receive one within your initial knee replacement operation.
- Crepitus.

- Damage to blood vessels, temporary or permanent nerve damage resulting in pain or numbness to the affected limb (nerve damage is more often a result of surgical trauma rather than being caused by the device).
- Sudden drop in blood pressure intra-operatively due to the use of bone cement.
- Valgus-varus deformity (position of joint (or limb) is angled outward or inward).
- Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarction, cerebrovascular event.
- Swelling/oedema, bleeding, haematoma (localised bleeding outside of blood vessel), hemarthrosis (bleeding into the joint cavity).
- Delayed wound healing, and other wound healing complications.
- Periarticular calcification (accumulation of calcium in the joint) or ossification (bone formation) with loss of joint mobility.
- Release of metal ions due to leaching or corrosion of metal components.
- Allergic reactions.
- The device has limited ability to withstand excessive activity levels and extra care should be taken during the post-operative recovery period. You should avoid putting excessive demands on the joint replacement such as excessive activity, failure to control your body weight and trauma, which have been implicated with early failure by loosening, fracture, and/or wear of the device. The device does not replace normal healthy bone, and the device can break or be damaged as a result of strenuous activity or trauma.
- Loosening of the components can result in increased production of wear particles, as well as accelerate damage to the bone, making successful revision surgery more difficult. Postoperative care is important and you are cautioned to govern activities accordingly, protecting the joint replacement from unreasonable stresses.

To help reduce the likelihood of residual risks, follow the instructions provided by your surgeon with respect to self-care, follow-up care, rehabilitation, physiotherapy and any other treatment in relation to your operation.

You should contact your healthcare professional if you believe that you are experiencing undesirable side-effects related to the device or its use, or if you are concerned about the risks. This document is not intended to replace a consultation with your surgeon, if needed.

## **Warnings, precautions or measures with regard to exposure to reasonably foreseeable external environments including risks of device interaction with other equipment**

The SAIPH® Knee System has been tested for interaction with Magnetic Resonance Imaging (MRI) equipment. General recommendations for MR parameters which should not expose you to additional hazard during MRI are provided in the Instructions for Use which your healthcare provider can access here: '<https://info.matortho.com/reg/400-480.pdf>.'

You should make your doctor aware of any other implants or known health considerations that may require reduction of the MRI limits. There are several different manufacturers and generations of MRI systems available and MatOrtho® cannot make any claims regarding the safety of MatOrtho® implants and devices with any specific MRI system.

## **Symptoms that could indicate that the device is malfunctioning**

Unexpected pain, stiffness, instability or a general worsening of symptoms around your knee especially when compared to your preoperative or earlier timepoint post-operative condition may indicate that there is a problem with your joint replacement. An inflamed, hot, itchy or red area around your knee or wound may indicate infection.

## **Precautions and measures that should be taken if the performance of the device changes or if any symptoms mentioned in this leaflet are experienced**

A general worsening of the performance of the replaced knee as compared to the preoperative condition may indicate that there is a problem with the joint replacement. Serious incidents that occurs in relation to the device should be reported to the appropriate country authority:

- Australia: Therapeutic Goods Administration (TGA). Information for this can be found on the Therapeutic Goods Administration's website: <https://www.tga.gov.au/medical-device-incident-reporting-investigation-scheme-iris>.
- UK: <https://yellowcard.mhra.gov.uk/>
- EU (Authorised Representative):
  - > [AR@matortho.com](mailto:AR@matortho.com)
  - > <https://www.hpra.ie/about-us/patient-forum/patient-forum-resources>
  - > <https://www.hpra.ie/safety-information/how-we-monitor-safety/medical-devices/reporting-safety-issues-for-devices/reporting-device-safety-issues-for-healthcare-professionals>

- Canada:  
<https://www.canada.ca/en/health-canada/services/drugs-health-products/reports-publications/medeffect-canada.html>

The healthcare professional should inform the manufacturer of any serious adverse events by emailing 'per.admin@matortho.com.'

### **Precautions and other measures that should be taken at, or near, the end of the expected device lifetime**

The precautions, risks, and warnings described above are appropriate throughout the lifetime of the device.

### **Other circumstances in which you should contact a health professional in relation to the operation of the device**

No further specific circumstances beyond those discussed in this leaflet.

### **Manufacturer**

For further information, please visit MatOrtho® patient information website:  
[www.matortho.com](http://www.matortho.com)

or contact MatOrtho® at:

MatOrtho Limited, 19/20 Mole Business Park, Randalls Road, Leatherhead, Surrey, KT22 7BA, United Kingdom.

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